









Model	Vismo series (PVM-47XX)	Vismo series (PVM-27XX)	Vital Sign Monitor series (SVM-7XXX)	SVM-75XX/76XXseries	OLG-3800 series	OLV-4201series
Appearance						
Chlorhexidine gluconate solution: 0.5%	✓	✓	✓	✓	✓	✓
Benzethonium chloride solution: 0.2%	✓	✓	✓	✓	✓	✓
Glutaraldehyde solution: 2.0%	✓	✓	✓	✓		
Benzalkonium chloride: 0.2%	✓	✓	✓	✓	✓	✓
Alkyldiaminoethylglycine hydrochloride: 0.5%	✓	✓	✓	✓	✓	✓
Phtharal: 0.55%	✓		✓	✓	✓	
Phenol: 1.56%	✓(1.8~2.3w/v)		✓(1.8~2.3w/v)	✓	✓(2.2%)	✓(2.2%)
Isopropyl alcohol: 70% by vol	✓		✓	✓	✓	✓
Ethanol: 76.9 to 81.4 vol% at15°C (59°F)	✓	✓	✓	✓	✓	✓
Neutral detergent (Cleaning)	✓	✓	✓	✓	✓	✓

sodium hypochlorite

Model	ZS-530P	ZS-611P	ZS-620P	ZS-630P	ZS-640P	GZ-130/140
Appearance						
Chlorhexidine gluconate solution: 0.5%	✓	✓	✓	✓	✓	✓
Benzethonium chloride solution: 0.2%	✓	✓	✓	✓	✓	✓
Glutaraldehyde solution: 2.0%	✓	✓	✓	✓	✓	✓
Benzalkonium chloride: 0.2%	✓	✓	✓	✓	✓	✓
Alkyldiaminoethylglycine hydrochloride: 0.5%	✓	✓	✓	✓	✓	✓
Phtharal: 0.55%		✓	✓	✓	✓	✓
Phenol: 1.56%		✓	✓	✓	✓	✓
Isopropyl alcohol: 70% by vol		✓	✓	✓	✓	✓
Ethanol: 76.9 to 81.4 vol% at 15°C (59°F)	✓	✓	✓	✓	✓	✓
Neutral detergent (Cleaning)	✓	✓	✓	✓	✓	✓

sodium hypochlorite