









Model	CO2 sensor TG-980P	CO2 sensor TG-920P	CO2 sensor TG-900P	SpO2 sensor TL-201T	SpO2 sensor TL-631T	SpO2 sensor TL-220T	SpO2 sensor TL-651T	Cuff inflation bag	Cuff cover
Appearance									
Chlorhexidine gluconate solution: 0.5%		✓✘	✓✘						
Benzethonium chloride solution: 0.2%		✓✘	✓✘		✓	✓			
Glutaraldehyde solution: 2.0%				✓	✓	✓	✓		
Benzalkonium chloride: 0.2%		✓✘	✓✘		✓	✓			
Alkyldiaminoethylglycine hydrochloride: 0.5%		✓✘	✓✘	✓	✓	✓			
Phtharal: 0.55%							✓		
Phenol: 1.56%									
Isopropyl alcohol: 70% by vol									
Ethanol: 76.9 to 81.4 vol% at15°C (59°F)	✓	✓	✓	✓	✓	✓	✓		
Neutral detergent	✓	✓	✓					✓	✓
sodium hypochlorite								✓	

✓: Applicable

✓✘: Applicable only for photo detector and light emitter