

Multi-link Cable

YS-096P5

General

The YS-096P5 Multi-link Cable connects the multi-link sockets on a QI-600P interface unit and a QI-631P, QI-632P, QI-671P or QI-672P interface installed in a BSM-6000A/K series bedside monitor.

Safety Information

For safety information, refer to the Operator's Manual or Installation Guide of the device using this product.

This Safety and Performance Information is an extract from the general and safety information sections of the most recent edition of Operator's Manual or Installation Guide. Therefore, the contents of your Operator's Manual or Installation Guide may differ from those of this Safety and Performance Information. For detailed operating procedures, follow the instructions of your Operator's Manual or Installation Guide.

 Manufacturer
NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku,
Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>

 European Representative
NIHON KOHDEN EUROPE GmbH
Raiffisenstrasse 10, D-61191 Rosbach, Germany
Phone +49 6003-827-0 Fax +49 6003-827-599