


MK-3CN HEMATOLOGY CONTROL+CRP 2.5mL


ASSAY SHEET

LOT	MK3C23YN
	2024-02-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	7.56	± 0.80
RBC	10 ⁶ /μL	4.68	± 0.20
HGB	g/dL	13.77	± 0.50
HCT	%	45.3	± 3.0
MCV	fL	96.9	± 7.0
MCH	pg	29.4	± 2.4
MCHC	g/dL	30.4	± 3.0
PLT	10 ³ /μL	279.8	± 50.0
MPV	fL	8.7	± 2.0
RDW-CV	%	16.2	± 3.0
LY%	%	34.15	± 10.00
MO%	%	2.51	± 5.00
GR%	%	63.34	± 10.00
LY	10 ³ /μL	2.58	± 1.00
MO	10 ³ /μL	0.19	± 0.40
GR	10 ³ /μL	4.79	± 1.00
CRP	Precision	≤ 10.0%	
		(16.30-30.30 mg/L)	

Manual Registration	
Lot	23Y
Level	Normal
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN

0694-901766D