


**MK-3CH** HEMATOLOGY CONTROL+CRP 2.5mL


**ASSAY SHEET**

<b>LOT</b>	<b>MK3C242H</b>
	<b>2024-05-05</b>
<b>OPEN</b>	<b>14DAYS</b>

**MEK-1303 Barcode Entry**

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 <sup>3</sup> /μL	20.07	± 2.20
RBC	10 <sup>6</sup> /μL	5.75	± 0.25
HGB	g/dL	18.97	± 0.70
HCT	%	60.8	± 4.0
MCV	fL	105.8	± 7.0
MCH	pg	33.0	± 2.8
MCHC	g/dL	31.2	± 3.0
PLT	10 <sup>3</sup> /μL	575.6	± 70.0
MPV	fL	8.5	± 3.0
RDW-CV	%	16.6	± 3.0
LY%	%	15.43	± 10.00
MO%	%	1.44	± 5.00
GR%	%	83.13	± 10.00
LY	10 <sup>3</sup> /μL	3.10	± 1.80
MO	10 <sup>3</sup> /μL	0.29	± 0.70
GR	10 <sup>3</sup> /μL	16.68	± 2.00
CRP	Precision	≤ 10.0%	
		(42.50-78.90 mg/L)	

Manual Registration	
Lot	<b>242</b>
Level	<b>High</b>
Control	<b>MK-3C</b>


<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**

 1-31-4 NISHIOCHIAI, SHINJUKU-KU  
 TOKYO 161-8560 JAPAN

**0694-901767E**