


MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS243
	2024-04-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



		MEK -6400 -6410 -6420	MEK -6500 -6510	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305	MEK -6450	MEK -6550
WBC	10 ³ /μL	9.1	9.1	9.3	9.2	9.0	9.43	9.19	9.1	9.1
RBC	10 ⁶ /μL	4.69	4.69	4.68	4.65	4.67	4.74	4.67	4.69	4.69
HGB	g/dL	13.9	14.0	13.7	13.7	13.6	13.97	13.87	13.9	14.0
HCT	%	40.5	40.5	42.2	41.0	39.8	47.2	46.4	40.5	40.5
MCV	fL	86.4	86.4	90.1	88.1	85.2	99.5	99.4	86.4	86.4
MCH	pg	29.6	29.9	29.3	29.5	29.1	29.5	29.7	29.6	29.9
MCHC	g/dL	34.3	34.6	32.5	33.4	34.2	29.6	29.9	34.3	34.6
PLT	10 ³ /μL	253	253	242	245	244	263.1	270.0	253	253
MPV	fL	7.4	9.1	7.4	7.9	7.4	9.6	9.5	7.4	9.1
RDW-CV *	%	15.7	15.7	14.6	15.9	14.9	15.5	15.5	15.7	15.7

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	169
FL	87
SD	109
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901780M